

DOLY

LIGHTWEIGHT COMPACT MANUAL WHEELCHAIR



LIGHTWEIGHT
COMPACT



REMOVABLE
SANDWICH CUSHION



MULTICOLOR
AVAILABLE



STREAMLINED
MODEL

DOLY

LIGHTWEIGHT COMPACT MANUAL WHEELCHAIR

ELEGANT AND STYLISH

Streamlined body, available in multiple colors, perfectly blending aesthetics and functionality.



REAR WHEELS

14 inch / **20** inch

MAX CAPACITY

120 kg

FOLDING WIDTH

25.5 cm (14") / **27** cm (20")

LENGTH×WIDTH×HEIGHT

82.5 cm · **59.5** cm · **90** cm (14")
93 cm · **64** cm · **90** cm (20")

G./N. WEIGHT

11.8 kg / **9.8** kg (14")
12.8 kg / **10.8** kg (20")



FOLDING DESIGN

Foldable backrest and collapsible frame for easy storage and portability

LIGHTWEIGHT FREEDOM

Featuring an ultra-light aluminum alloy frame, weighs around 10kg, move freely and lightly



DUAL BRAKES

Safety is paramount with Doly's dual braking system, includes handbrakes and cable brakes for reliable stopping



OPTIONAL SIZES

Adapt to your environment with interchangeable rear wheels and convenient handbrakes for easy adjustment

