

LIFT




LIFT Y401

STEEL

/ FRAME	Steel, powder coating
/ CHASSIS LENGTH	1070mm
/ CHASSIS WIDTH	600mm - 800mm
/ LIFTING HEIGHT	830mm - 1850mm
/ HANDLE HEIGHT	1000mm
/ CASTOR HEIGHT	90mm
/ MOTOR	DC 24V 8000N
/ BATTERY	DC 24V 5Ah, rechargeable battery Approx. 60 lifts per full charge
/ CHARGER	24V-1A
/ CHARGING TIME	6 hours per time
/ FRONT CASTOR	3" solid tire, plastic wheel optional
/ REAR CASTOR	4" solid tire with locking, plastic wheel optional
/ MAX. USER WEIGHT	150KGS (23.62 stone)
/ ALTERNATIVES	Purpose designed sling
/ PACKAGE	1260×630×350mm
/ G./N. W.	51 / 43KGS



Y402 STANDING



Y403 FOLDED HOME-USE



Y404 MANUAL

ITEM NO.	Chassis Length (mm)	Chassis Width (mm)	Lifting Height (mm)	Handle Height (mm)	Castor Height (mm)	Motor	Charger	Front Caster Spec.	Rear Caster Spec.	Max. User Weight (kg)	G./N. Weight (kg)	Carton Size (mm)
Y402	1070	600-800	1295-1960	1000	90	DC 24V 8000N	24V-1A	3"	4"	150	56.5/48.5	1240x630x350
Y403	940	535-840	770-1830	1080	90	DC 24V 8000N	24V-1A	3"	4"	150	47.0/39.5	1260x570x350
Y404	1160	600-800	910-1900	1000	90	DC 24V 8000N	24V-1A	3"	4"	150	51.0/43.0	1240x630x350

LIFT Y407 STEEL

/ FRAME	Steel, powder coating
/ CHASSIS LENGTH	1200mm
/ CHASSIS WIDTH	600mm - 1065mm
/ LIFTING HEIGHT	840mm - 1870mm
/ HANDLE HEIGHT	1170mm
/ CASTOR HEIGHT	115mm
/ MOTOR	DC 24V 8000N
/ BATTERY	DC 24V 2.2AH, interchangeable feeding battery Approx. 60 lifts per full charge
/ CHARGER	24V-1A
/ CHARGING TIME	5 hours per time
/ FRONT CASTOR	3" solid tires
/ REAR CASTOR	4" solid tires with locking
/ PUSH HANDLE	Push handle with rounded corners
/ MAX. USER WEIGHT	150KGS (23,62 stone)
/ ALTERNATIVES	Purpose designed sling(3 kinds) Standard Mode / Stand Aid Mode Interchangeable
/ PACKAGE	1240×630×360mm
/ G./N. W.	48/40KGS



The hand lever controls the base to open and close



The double rear wheels are equipped with safety lock device



The hanger can rotate 360 ° at will

LIFT Y410

STEEL

NEW

/ FRAME	Steel, powder coating
/ SEAT TO FLOOR HEIGHT	420mm-570mm
/ TOTAL LENGTH	730mm
/ TOTAL WIDTH	560mm
/ TOTAL HEIGHT	810mm
/ FOLDED / DISMANTLED LENGTH	360mm
/ FOLDED / DISMANTLED WIDTH	560mm
/ FOLDED / DISMANTLED HEIGHT	680mm
/ LEG HEIGHT / FRAME HEIGHT	900mm
/ SUPPLY VOLTAGE	24V
/ MAXIMUM LOAD	120kg
/ LIFTING RANGE	150mm
/ BATTERY CAPACITY	24V4.4Ah
/ TOTAL WEIGHT	35kg
/ PACKAGE	600×430×760mm



Effortlessly and easily folds down into three sections



One key lifting, easy operation



Metal safety buckle



Silent universal wheel with brake



Steel pedals, strong and durable