


ROLLATOR














SPiRiT X3

ALUMINUM

/ FRAME	Aluminum
/ WHEELS	9", PU tires/ 8", TPR tires
/ SEAT	Supplied
/ SHOPPING BASKET	Supplied
/ HANDBRAKES	With locking mechanism
/ COLOR	Multiple colors can be customized 
/ PACKAGE	228pcs/ 20" GP, 591pcs/ 40" HQ


NEW



ITEM NO.	Seat Height (cm)	Width Inside Hand Grip (cm)	Overall Width (cm)	Overall Length (cm)	Seat Width (cm)	Seat Depth (cm)	Front Caster Spec.	Height Adjustable (cm)	Weight Capacity (kg)	G./N. Weight (kg)	Carton Size (cm)
											
SPiRiT X3	51	47.5	60.5	72	/	/	9"	68.5-78.7	130	8.2/7.7	73x25.55x77
	55	47.5	60.5	72	/	/	9"	68.2-82.8	130	8.4/7.9	73x25.5x79.5
	58	47.5	60.5	72	/	/	9"	73.7-91.4	130	8.5/8.0	73x25.5x82.5
	64	47.5	60.5	72	/	/	9"	76.8-97.2	130	8.7/8.2	73.5x25.4x87.5












SPiRiT X4

CARBON FIBER

/ FRAME	Carbon Fiber
/ WHEELS	8", TPR tires
/ SEAT	Supplied
/ SHOPPING BASKET	Supplied
/ HANDBRAKES	With locking mechanism
/ COLOR	Multiple colors can be customized 
/ PACKAGE	228pcs/ 20" GP, 591pcs/ 40" HQ

NEW



ITEM NO.	Seat Height (cm)	Width Inside Hand Grip (cm)	Overall Width (cm)	Overall Length (cm)	Seat Width (cm)	Seat Depth (cm)	Front Caster Spec.	Height Adjustable (cm)	Weight Capacity (kg)	G./N. Weight (kg)	Carton Size (cm)
											
SPiRiT X4	51	47.5	60.5	72	/	/	9"	70-78.7	130	7.6/7.1	73x25.5x77
	55	47.5	60.5	72	/	/	8"	69.7-82.8	130	7.7/7.2	73x25.5x79.5
	58	47.5	60.5	72	/	/	8"	75.2-91.4	130	7.8/7.3	73x25.5x82.5
	64	47.5	60.5	72	/	/	8"	76.8-97.2	130	7.9/7.4	73.5x25.4x87.5

YK7010

ALUMINUM

/ FRAME	Adjustable height
/ WHEELS	6"x1", TPR tires
/ SEAT	Padded with back support
/ SHOPPING BASKET	Supplied
/ HANDBRAKES	With locking mechanism
/ COLOR	Available in several colors
/ PACKAGE	216pcs/ 20" GP, 464pcs/ 40" GP, 520pcs/ 40" HQ



ITEM NO.	Seat Height (cm)	Width Inside Hand Grip (cm)	Overall Width (cm)	Overall Length (cm)	Seat Width (cm)	Seat Depth (cm)	Front Caster Spec.	Height Adjustable (cm)	Weight Capacity (kg)	G./N. Weight (kg)	Carton Size (cm)
YK7010	46	62	59	45	34	56	6"	85-96	110	8.6/6.6	62.5x23x85.5

YK7030

ALUMINUM

/ FRAME	Adjustable height
/ WHEELS	8"x1"
/ SEAT	ABS
/ SHOPPING BASKET	Supplied
/ HANDBRAKES	With locking mechanism
/ COLOR	Available in several colors
/ PACKAGE	230pcs/ 20" GP, 485pcs/ 40" GP, 540pcs/ 40" HQ



ITEM NO.	Seat Height (cm)	Width Inside Hand Grip (cm)	Overall Width (cm)	Overall Length (cm)	Seat Width (cm)	Seat Depth (cm)	Front Caster Spec.	Height Adjustable (cm)	Weight Capacity (kg)	G./N. Weight (kg)	Carton Size (cm)
YK7030	46	60	67	45	34	52	8"	81-92	110	9.6/7.6	59x26x71

YK7050

STEEL

/ FRAME	Adjustable height
/ WHEELS	8"x2", TPR tires
/ SEAT	Plastic with tray
/ SHOPPING BASKET	Supplied
/ HANDBRAKES	With locking mechanism
/ COLOR	Available in several colors
/ PACKAGE	260pcs/ 20" GP, 540pcs/ 40" GP, 640pcs/ 40" HQ



ITEM NO.	Seat Height (cm)	Width Inside Hand Grip (cm)	Overall Width (cm)	Overall Length (cm)	Seat Width (cm)	Seat Depth (cm)	Front Caster Spec.	Height Adjustable (cm)	Weight Capacity (kg)	G./N. Weight (kg)	Carton Size (cm)
YK7050	61	43	65	65	43	/	8"	83-96	110	15.0/13.0	59x26x71

YK7060

STEEL

/ FRAME	Adjustable height
/ WHEELS	8"x2", TPR tires
/ SEAT	/
/ SHOPPING BASKET	Supplied
/ HANDBRAKES	With locking mechanism
/ COLOR	Available in several colors
/ PACKAGE	336pcs/ 20" GP, 684pcs/ 40" GP, 756pcs/ 40" HQ (2pcs/carton)



ITEM NO.	Seat Height (cm)	Width Inside Hand Grip (cm)	Overall Width (cm)	Overall Length (cm)	Seat Width (cm)	Seat Depth (cm)	Front Caster Spec.	Height Adjustable (cm)	Weight Capacity (kg)	G./N. Weight (kg)	Carton Size (cm)
YK7060	/	55	66	/	/	/	8"	82-96	/	9.1/7.1	66x36x75