

GEMINI

HIGH ENDURANCE, HEAVY-DUTY LARGE SCOOTER



HIGH-CAPACITY
LEAD-ACID BATTERY



POWERFUL
MOTOR



FRONT AND REAR
SUSPENSION



PERSONALIZED
SEATING SYSTEM

GEMINI

BRAND NEW PRODUCT LAUNCH

JOURNEY OF FREEDOM

Ultra-long range, superior weight capacity, excellent shock absorption, adjustable speed, and a journey of your choice.



MAX SPEED

15 km/h

TRAVEL MILEAGE

50 km

MAX SLOPE

10°

MAX CAPACITY


150 kg

LENGTH×WIDTH×HEIGHT

136.5 cm · **62.5** cm · **133** cm

TOTAL MASS(BATTERY INCLUDED/EXCLUDED)

135 kg / **100** kg




HIGH-CAPACITY BATTERY

24V 75Ah stable and durable lead-acid battery, travel mileage up to 50 km.

POWERFUL MOTOR

450W/600W motor: speed up to 15 km/h, unfazed by slopes.



TWO CUP HOLDERS

The handlebar housing includes two built-in cup holders for convenient storage.

CUSTOM SEAT SYSTEM

360° rotation adjustable seat; flip-back armrests; adjustable headrest.



SUSPENSION

13" pneumatic tires with aluminum alloy wheels for superior obstacle-crossing; front and rear shock absorbers for enhanced comfort.

