

# SCOOTER

---



# VIRGO-D VIRGO-F

DETACHABLE VERSION

NEW



/ CONTROLLER	Dynamic 50A
/ MOTOR	270W, electromagnetic brake
/ SLOPE	6°
/ MAX SPEED	8 KM/H
/ TRAVEL MILEAGE	16 KM
/ BATTERY	Lead acid 24V 12Ah
/ WHEEL	8"x2.5", black PU tire
/ CHARGER	24V 2A
/ SEAT	PVC with sponge
/ COLOR	Champagne/ Red/ Blue
/ PACKAGE	227pcs/40"HQ
/ SUSPENSION	VIRGO-D: / ; VIRGO-F:Rear




Delta tiller with dashboard



360 ° rotatable ergonomic seat



Equipped with storage basket



Rearview mirror optional



VIRGO-D  
Easily disassembled into five pieces



VIRGO-F  
Easily disassembled into four pieces

ITEM NO.	Total Width (cm)	Total Height (cm)	Total Length (cm)	Seat Height (cm)	Seat Width (cm)	Seat Depth (cm)	G./N. Weight (kg)	Weight Capacity (kg)	Carton Size (cm)
VIRGO-D	50	85.5-89	103	49,5-53	44	39	49.0/45.0	136	113x50,5x53

# LIBRA-D

/CONTROLLER	Dynamic 50A
/MOTOR	270W, electromagnetic brake
/SLOPE	6°
/MAX SPEED	8 KM/H
/TRAVEL MILEAGE	16 KM
/BATTERY	Lead acid 24V 15Ah
/WHEEL	8"×2.5", black PU tire
/CHARGER	24V 2A
/SEAT	PVC with sponge
/COLOR	Champagne/ Red/ Blue
/PACKAGE	227pcs/40"HQ



Delta tiller with dashboard




360 ° rotatable ergonomic seat



Easily disassembled into five pieces



Equipped with storage basket



Rearview mirror optional



Position light

ITEM NO.	Total Width (cm)	Total Height (cm)	Total Length (cm)	Seat Height (cm)	Seat Width (cm)	Seat Depth (cm)	G./N. Weight (kg)	Weight Capacity (kg)	Carton Size (cm)
LIBRA-D	48.5	88	108.3	51-54.5	43.5	38	39.5/35.5	136	113x50.5x53



# LIBRA-DX

NEW

/CONTROLLER	Dynamic 50A
/MOTOR	270W, electromagnetic brake
/SLOPE	6°
/MAX SPEED	8 KM/H
/TRAVEL MILEAGE	25 KM
/BATTERY	Lead acid 24V 25Ah
/WHEEL	8"x2,5", black PU tire
/CHARGER	24V 2A
/SEAT	PVC with soft sponge
/COLOR	Champagne/ Red/ Blue
/PACKAGE	227pcs/40"HQ



Delta tiller with dashboard



360 ° rotatable ergonomic seat




Easily disassembled into five pieces



Equipped with storage basket



Position light



Front and rear suspension


ITEM NO.	Total Width (cm)	Total Height (cm)	Total Length (cm)	Seat Height (cm)	Seat Width (cm)	Seat Depth (cm)	G./N. Weight (kg)	Weight Capacity (kg)	Carton Size (cm)
LIBRA-DX	48.5	88	108.3	53.5-57	44	41	57.1/52.9	136	113x50.5x53

# SAGITTA-D

/ CONTROLLER	Dynamic DR90
/ MOTOR	320W, electromagnetic brake
/ SLOPE	9°
/ MAX SPEED	10 KM/H
/ TRAVEL MILEAGE	38 KM/ 50KM
/ BATTERY	24V 38Ah/ 24V 52Ah
/ WHEEL	3.00-4, air tyre
/ CHARGER	24V 4A/ 24V 6A
/ SEAT	PVC Leather with sponge
/ PACKAGE	36pcs/20"GP, 78pcs/40"GP



Intuitive control panel



Comfort adjustable seating



Easily disassembled into five pieces



Equipped with storage basket



Brake light/indicators



Headlight

ITEM NO.		Total Width (cm)	Total Height (cm)	Total Length (cm)	Seat Height (cm)	Seat Width (cm)	Seat Depth (cm)	G./N. Weight (kg)	Weight Capacity (kg)	Carton Size (cm)
SAGITTA-D	38Ah	58.3	96.5	125.5	50-54	45	45.5	77.6/89.6	150	Body: 133x63x63
	52Ah	58.3	96.5	125.5	50-54	45	45.5	83.0/93.6	150	Seat: 60x55x39