

BUGGY



YK8001

ALUMINUM

/ FRAME	Light weight and luxury design
/ FRONT WHEEL	8"x2", PU tire
/ REAR WHEEL	12", PU tire
/ PUSHBAR	Height adjustable
/ FOOTREST	Adjustable angle and height footrest
/ SEAT AND BACKREST	Changeable seat width and depth; Comfortable soft and adjustable headrest
/ FRONT FORK	Plastic
/ SEAT SYSTEM	Removable, two-way adjustable and reclining
/ PACKAGE	121pcs/ 20" GP, 230pcs/ 40" HQ
/ SPECIAL	Foldable chair frame for easy transportation and storage; Mesh bag for small belongings



ITEM NO.	SIZE	Overall Length (cm)	Overall Width (cm)	Overall Height (cm)	Overall Width Folded	Seat Width (cm)	Seat Depth (cm)	Back Height (cm)	Floor to Seat Height	Weight Capacity (kg)	Front Caster Spec.	Rear Wheel Spec.	G./N. Weight (kg)	Carton Size (cm)
YK8001		83	64	116	/	28-32	30	/	51.5	60	8"	12"	22.4/19.6	82x66x45
		83	64	116	/	32-40	30	/	51.5	60	8"	12"	22.4/19.6	82x74x45

YK9140

ALUMINUM

/ FRAME	Single cross brace, powder coating
/ FRONT WHEEL	6", PU tire
/ REAR WHEEL	16", PU tire
/ ARMREST	Height adjustable armrest, PU armpad
/ FOOTREST	Elevating legrest
/ SEAT AND BACKREST	Breathable mesh fabric
/ FRONT FORK	Aluminum
/ FOOTPLATE	Plastic
/ PACKAGE	100pcs/ 20" GP, 200pcs/ 40" HQ
/ OPTION	Attendant brakes; Tray table
/ SPECIAL	Backrest reclining; Seat tilting



ITEM NO.	SIZE	Overall Length (cm)	Overall Width (cm)	Overall Height (cm)	Overall Width Folded	Seat Width (cm)	Seat Depth (cm)	Back Height (cm)	Floor to Seat Height	Weight Capacity (kg)	Front Caster Spec.	Rear Wheel Spec.	G./N. Weight (kg)	Carton Size (cm)
YK9140	14"	98	48	100	33	36	38	46	42	≤90	6"	16"	26.0/24.0	80×40×64
	16"	98	53	100	33	40	38	46	42	≤90	6"	16"	27.0/25.0	80×40×68
	18"	112	60	110	33	46	40	54	45	≤90	6"	16"	30.0/27.0	95×34×95