

COMMODE

CHAIR



YK4010 YK4030

STEEL

/ FRAME	Powder coating
/ REAR WHEEL	With brake
/ ARMREST	Full length padded armrest
/ FOOTREST	Swing away footrest
/ SEAT AND BACKREST	Removable padded seat and PVC backrest
/ COLOR	Plastic



YK4010



YK4030

ITEM NO.	Seat Width (cm)	Overall Height (cm)	Seat Height (cm)	Front Caster Spec.	Rear Wheel Spec.	Weight Capacity (kg)	G./N. Weight (kg)	Carton Size (cm)
YK4010	44	92	52	5"	5"	115	15.5/13.0	59x56x60
YK4030	44	95	52	5"	5"	115	14.5/12.5	58x52.5x54.5

YK4120

STEEL

/ FRAME	Powder coating
/ ARMREST	HDPE + S1010
/ TIPS	Rubber+A3
/ SEAT COVER	PP
/ BUCKET SET	PP
/ COLOR	Grey
/ CAPACITY	114kg/ 250lbs



ITEM NO.	Seat Width (cm)	Overall Height (cm)	Seat Height (cm)	Front Caster Spec.	Rear Wheel Spec.	Weight Capacity (kg)	G./N. Weight (kg)	Carton Size (cm)
YK4120	53	63.5-73.5	42-52	/	/	130	12.0/11.0	62x35x50

YK4020

STEEL

/ FRAME	Single cross brace, powder coating
/ FRONT WHEEL	8"x2", PU tires
/ REAR WHEEL	24" spoke wheels, PU tire
/ ARMREST	Detachable, full length armrest, PVC armpad
/ FOOTREST	Swing away footrest
/ SEAT AND BACKREST	U shape seat, PVC backrest
/ COMMODE	Plastic
/ FRONT FORK	Steel
/ FOOT PLATE	Plastic
/ PACKAGE	98pcs/ 20" GP, 228pcs/ 40" HQ



ITEM NO.	Seat Width (cm)	Overall Height (cm)	Seat Height (cm)	Front Caster Spec.	Rear Wheel Spec.	Weight Capacity (kg)	G./N. Weight (kg)	Carton Size (cm)
YK4020	45	92	53.5	8"	24"	115	24.1/21.4	82x35x94

YK4050

STEEL

/ FRAME	Single cross brace, powder coating
/ FRONT WHEEL	8"x2", PU tires
/ REAR WHEEL	24" spoke wheels, Pneumatic / PU tires
/ ARMREST	Flip back with desk length
/ FOOTREST	Swing away footrest (Detachable)
/ SEAT AND BACKREST	Removable padded seat and Nylon backrest
/ COMMODE	Plastic
/ FRONT FORK	Steel
/ FOOT PLATE	Plastic
/ PACKAGE	135pcs/ 20" GP, 270pcs/ 40" GP, 320pcs/ 40" HQ



ITEM NO.	Seat Width (cm)	Overall Height (cm)	Seat Height (cm)	Front Caster Spec.	Rear Wheel Spec.	Weight Capacity (kg)	G./N. Weight (kg)	Carton Size (cm)
YK4050	46	91	54	8"	24"	120	23.0/21.0	82x35x94
	51	91	54	8"	24"	120	24.5/22.5	82x35x94