

# POWER

LIFT CHAIR




# VITTO-1M

/STYLE	Classic
/MOTOR	Single Motor
/FRAME	Constructed with a mixed hardwood frame
/UPHOLSTERY	Tissue / PU cover
/SEAT	High density foam
/FEATURES	Electric reclining and lifting mechanism ; One remote for all the operations; Easily assembled in less than 15 minutes;
/COLOR OPTIONS	Beige / Gray / Black / Optional
/TOTAL WIDTH	800 mm
/SEAT WIDTH	540 mm
/ARMREST HEIGHT	140 mm
/SEAT HEIGHT	505 mm
/TOTAL HEIGHT	995 mm
/TOTAL LENGTH	1564 mm
/SEAT DEPTH	425 mm
/BACKREST HEIGHT	780 mm
/G.W.	52 kgs
/WEIGHT CAPACITY	150 kg
/WARRANTY	1 Year
/OPTIONAL	Dual motor; Massage; Heat
/PACKAGE	42pcs/20"GP,120pcs/40"HQ



TISSUE 221-1



TISSUE 221-14



TISSUE 221-23



BL0720-34



BL0720-37



YKM-03



YKPU-01



YKPU-02




YKPU-03

# YK3001-1M

/STYLE	Intelligence
/MOTOR	Single Motor
/FRAME	Reinforced and thickened steel pipes
/UPHOLSTERY	Tissue / PU cover
/SEAT	Filed high-density rebound sponge
/FEATURES	Electric reclining and lifting mechanism ; One remote for all the operations; Easily assembled in less than 15 minutes;
/COLOR OPTIONS	Beige / Gray / Black / Optional
/TOTAL WIDTH	800 mm
/SEAT WIDTH	550 mm
/ARMREST HEIGHT	125 mm
/SEAT HEIGHT	520 mm
/TOTAL HEIGHT	995 mm
/TOTAL LENGTH	1637 mm
/SEAT DEPTH	460 mm
/BACKREST HEIGHT	760 mm
/G.W.	45 kgs
/WEIGHT CAPACITY	150 kg
/WARRANTY	1 Year
/OPTIONAL	Dual motor; Massage; Heat
/PACKAGE	70pcs/20"GP, 177pcs/40"HQ



TISSUE 221-1



TISSUE 221-14



TISSUE 221-23



BL0720-34



BL0720-37



YKM-03



YKPU-01



YKPU-02



YKPU-03