

SHOWER

CHAIR



YK3010

/ FRAME	Height adjustable
/ TUBE	Aluminum
/ SEAT AND BACKREST	High-grade PVC; U-shaped cut



ITEM NO.	Seat Width (cm)	Seat Depth (cm)	Back Height (cm)	Floor to Seat Height (cm)	G./N. Weight (1pcs) (kg)	G./N. Weight (2pcs) (kg)	Carton Size 1 pcs per carton (cm)	Carton Size 2 pcs per carton (cm)
YK3010	50	30	/	43	2.4/1.9	5.0/4.0	60×11.5×36	60×18×36

YK3020

/ FRAME	Height adjustable with backrest
/ TUBE	Aluminum
/ SEAT AND BACKREST	High-grade PVC; U-shaped cut



ITEM NO.	Seat Width (cm)	Seat Depth (cm)	Back Height (cm)	Floor to Seat Height (cm)	G./N. Weight (1pcs) (kg)	G./N. Weight (2pcs) (kg)	Carton Size 1 pcs per carton (cm)	Carton Size 2 pcs per carton (cm)
YK3020	50	30	77	43	3.2/2.7	6.0/5.4	60×16×36	60×23×36