


CUSHION




L01 L09


/ MATERIAL	L01:High rebound material L09:Slow rebound material
/ WEIGHT CAPACITY	130 kg
/ DEPTH	430 mm
/ WIDTH	460 mm
/ SEAT THICKNESS	80 mm
/ LOWER THIGH SUPPORT THICKNESS	70 mm
/ INNER THIGH SUPPORT THICKNESS	45 mm
/ EXTERNAL THIGH SUPPORT HEIGHT	45 mm
/ PELVIC PROFILE HEIGHT	35 mm



Ergonomic design



Soft and comfortable fabric
Curve wrap around the hip



Double-sided honeycomb design



L01:High rebound material



L09:Slow rebound material



Waterproof hidden zipper

